







SUREBIOTM COVID-19 IgM/IgG

SV BIOTECH PVT. LTD.

CONCEPT TO COMMERCIALISATION BIOPHARMA VALUE CHAIN




TEST PRINCIPLE



IgM/IgG Rapid Test

- It is suited for POC Testing
- Provide information on previous infection
- Detects immune response & provide identification of likely immunity.
- Aids in verifying if people can be released from quarantine &/or return to work

Coronavirus is a single-stranded positive-sense RNA virus with an envelope of about 80 to 120 nm in diameter. Its genetic material is the largest of all RNA viruses and is an important pathogen of many domestic animals, pets, and human diseases. It can cause a variety of acute and chronic diseases. Common signs of a person infected with a coronavirus include respiratory symptoms, fever, cough, shortness of breath, and dyspnoea. In more severe cases, infection can cause pneumonia, severe acute respiratory syndrome, kidney failure, and even death. The 2019 new coronavirus, or "COVID-19", was discovered due to Wuhan Viral Pneumonia cases in 2019 and was named by the World Health Organization on January 12, 2020. WHO confirmed that COVID-19 can cause colds, the Middle East Respiratory Syndrome (MERS) and more serious diseases such as severe acute respiratory syndrome (SARS). This kit is helpful for the auxiliary diagnosis of coronavirus infection. The test results are for clinical reference only and cannot be used as a basis for confirming or excluding cases alone.



SV Biotech Offers COVID-19 IgG/IgM Antibody/Antigen RT-Kit from our collaborative/marketing Partner:

1. SureBio Covid-19 IgG/IgM (product launch within 2 month)
 2. P4Detect, South Korea (ISO 13485:2016/CE/ EU GMP)
 3. ImmunoQuick Covid-19 IgG/IgM (ICMR /CE/ISO 9001:2015/ISO 13485:2016).
- ** ImmunoScience India is Subsidiary of Immunoscience Inc

PRODUCT FEATURES

- Easier:** No special equipment needed Intuitive visual interpretation.
- Rapid:** Quick sampling by fingertip blood, result in 15 minutes.
- Accurate:** Qualitative lateral flow Immunoassay for the detection of Corona Virus IgM/IgG nucleoproteins in serum, plasma whole blood specimen.
- Application:** For suspicious patients with symptoms, mild symptoms, or even without symptoms, also for testing people with close contact of infected patients and people under quarantine control.

KIT STORAGE AND STABILITY

Store the kit at room temperature, 2-30°C / 36-86°F, out of direct sunlight. Kit materials are stable until the expiration date printed on the outer box. Do not freeze the kit.

SYMPTOMS



SV BIOTECH PVT. LTD. (OPC)

Flat No-801, Roystonia,
Bldg. K, Magarpatta City, Hadapsar,
PUNE-411028
Tel: +917887884108, 020-26911273/74
Email: Info@svbiotech.in
Web: www.svbiotech.in

